

Entire PD Channel Map

*Modified SSP

TOP of Cryo

USDaS					
PD Module	HB	SSP	SSPch	DAQch	
002-0047-FL34	Hamamatsu	SSP503	0-3	216	219
002-0008-IU54	Hamamatsu	SSP503	4-7	220	223
002-0058-FL24	Hamamatsu	SSP503	8-11	224	227
002-0063-IU19	Hamamatsu	SSP504	0-3	228	231
003-0026-FL07*	SensL-C1	SSP501	0-3	192	195
002-0014-IU26	Hamamatsu	SSP504	4-7	232	235
003-0024-FL33	SensL-C1	SSP501	4-7	196	199
003-0004-IU48	SensL-C1	SSP501	8-11	200	203
002-0041-FL36	Hamamatsu	SSP504	8-11	236	239
002-0036-IU47	SensL-C1	SSP502	0-3	204	207

MSDaS					
PD Module	HB	SSP	SSPch	DAQch	
002-0002-FL22	Hamamatsu	SSP601	0-3	240	243
002-0054-IU22	Hamamatsu	SSP601	4-7	244	247
002-0059-FL08	Hamamatsu	SSP601	8-11	248	251
002-0020-IU09	Hamamatsu	SSP602	0-3	252	255
002-0060-FL39	Hamamatsu	SSP602	4-7	256	259
ARAPUCA-2	Hamamatsu	SSP603	0-3	264	267
		SSP603	4-7	268	271
		SSP603	8-11	272	275
002-0055-FL40	Hamamatsu	SSP602	8-11	260	263
002-0013-IU01	Hamamatsu	SSP604	0-3	276	279
002-0011-FL15	Hamamatsu	SSP604	4-7	280	283
002-0031-IU02	Hamamatsu	SSP604	8-11	284	287

DSDaS					
PD Module	HB	SSP	SSPch	DAQch	
001-0003-FL01	SensL-C1	SSP401	0-3	144	147
002-0044-IU50	SensL-C1	SSP401	4-7	148	151
002-0039-FL29	SensL-A1	SSP401	8-11	152	155
003-0002-IU27	SensL-C1	SSP402	0-3	156	159
002-0025-FL25	SensL-C1	SSP402	4-7	160	163
003-0011-IU37	SensL-C1	SSP402	8-11	164	167
003-0048-FL42	SensL-C1	SSP403	0-3	168	171
002-0023-IU53	SensL-C1	SSP403	4-7	172	175
002-0038-IU35	SensL-C1	SSP403	8-11	176	179
002-0040-FLP06*	SensL-C1	SSP404	0-3	180	183

Bottom of Cryo

TOP of Cryo


USRaS					
PD Module	HB	SSP	SSPch	DAQch	
003-0031-IU20	SensL-A1	SSP301	0-3	96	99
002-0055-FL03	SensL-A1	SSP301	4-7	100	103
002-0020-IU31	SensL-A1	SSP301	8-11	104	107
ARAPUCA-1	Hamamatsu	SSP304	0-3	132	135
		SSP304	4-7	136	139
		SSP304	8-11	140	143
002-0042-IU52	SensL-A1	SSP302	0-3	108	111
002-0056-FL30	SensL-A1	SSP302	4-7	112	115
002-0047-IU17	SensL-A1	SSP302	8-11	116	119
002-0054-FL38	SensL-A1	SSP303	0-3	120	123
001-0039-IU51	SensL-A1	SSP303	4-7	124	127
003-0015-FL04	SensL-C1	SSP303	8-11	128	131

MSRaS					
PD Module	HB	SSP	SSPch	DAQch	
002-0049-IU16	SensL-A1	SSP201	0-3	48	51
001-0054-FL18	SensL-A1	SSP201	4-7	52	55
002-0035-IU13	SensL-A1	SSP201	8-11	56	59
002-0006-FL14	SensL-A1	SSP202	0-3	60	63
001-0044-IU18	SensL-A1	SSP202	4-7	64	67
002-0012-FL19	SensL-A1	SSP202	8-11	68	71
002-0027-IU12	SensL-A1	SSP203	0-3	72	75
002-0015-FL21	SensL-A1	SSP203	4-7	76	79
001-0052-IU14	SensL-A1	SSP203	8-11	80	83
003-0025-FL06	SensL-A1	SSP204	0-3	84	87

DSRaS					
PD Module	HB	SSP	SSPch	DAQch	
403-003-0063-IU28	SensL-A1	SSP101	0-3	0	3
403-003-0041-FL9	SensL-A1	SSP101	4-7	4	7
403-002-0001-IU15	SensL-A1	SSP101	8-11	8	11
403-003-0054-FLP12	SensL-A1	SSP102	0-3	12	15
403-001-0006-IU49	SensL-A1	SSP102	4-7	16	19
403-003-0064-FLP13	SensL-A1	SSP102	8-11	20	23
403-001-0061-IU04	SensL-A1	SSP103	0-3	24	27
403-001-0042-FLP4	SensL-A1	SSP103	4-7	28	31
403-001-0025-IU21	SensL-A1	SSP103	8-11	32	35
403-003-0020-FL5	SensL-A1	SSP104	0-3	36	39

Bottom of Cryo

SSP_Serial#



SSP_IP#

